

# Hockey Proff

2017



Hockey Proff is CE-marked and TÜV-tested

Hockey Proff is a Norwegian produced hockey sledge of high quality. The frame is produced from super light weight high-grade steel, with a specially manufactured anodized iceskate bracket in aluminum.



Hockey Proff with Nordic Lightning carbon sticks

The low seat height combined with a vacuum cast plastic seat with padding, makes you sit safe and stable on the sledge.

The entire sledge weighs only 5 kg.

With a main frame in 3 different sizes and seat in 3 different sizes, the Hockey Proff will fit everyone from small children to big adults.

The seat and iceskate bracket can be adjusted independently forwards and backwards to find the correct balance point on the hockey sledge. If more support is needed, it is possible to get an anatomical sportseat, which has a higher backrest, which also provides good side support.

In order to capture the different level of the players, Hockey Proff is delivered with skate brackets in 4 different widths; 20 mm, 34 mm, 50 mm and 90 mm. The better you become, the narrower skate bracket you can use.

### Standard equipment:

- Seat belts with fixlock buckle
- Footrest with fixlock buckle
- Iceskates including bracket



*Pictures shown may vary from actual products.*

## Technical data:

Model:	Short	Standard	Long
Art.no.	500045-K	500045	500045-L
Seat width	24-30cm	28-36cm	34-44cm
Seat depth	42cm	45cm	46cm
Back height	18cm	20cm	22cm
Leg length	63-96cm	93-117cm	93-132cm
Seat height	12,5cm	12,5cm	12,5cm
Max user weight	110kg	110kg	110kg
Total weight	4.8kg	5.6kg	5.8kg

## Accessories:

Art.no.	Description
500091-5	Bracket 5cm
503269	M2 seat belts with jacking mechanism, S, set
503270	M2 seat belts with jacking mechanism, M, set
503271	M2 seat belts with jacking mechanism, L, set
503168	Support Skates incl bracket, set
503141	Pushbar with mounting bracket, angle adjustable, removable
503158	Bracket for sports seat
500110	Sport Seat SW26, SD29, BH38 including padding
500111	Sport Seat SW30, SD35, BH44 including padding
500112	Sport Seat SW34, SD41, BH50 including padding
500113	Sport Seat SW38, SD47, BH56 including padding
500161	Hockey Stick Carbon Nordic Lightning w picks, L 70 cm, right
500162	Hockey Stick Carbon Nordic Lightning w picks, L 70 cm, left
500163	Hockey Stick Carbon Nordic Lightning w picks, L 72 cm, right
500164	Hockey Stick Carbon Nordic Lightning w picks, L 72 cm, left
500165	Hockey Stick Carbon Nordic Lightning w picks, L 74 cm, right
500166	Hockey Stick Carbon Nordic Lightning w picks, L 74 cm, left
500167	Hockey Stick Carbon Nordic Lightning w picks, L 76 cm, right
500168	Hockey Stick Carbon Nordic Lightning w picks, L 76 cm, left
500169	Hockey Stick Carbon Nordic Lightning w picks, L 78 cm, right
500170	Hockey Stick Carbon Nordic Lightning w picks, L 78 cm, left
500171	Hockey Stick Carbon Nordic Lightning w picks, L 80 cm, right
500172	Hockey Stick Carbon Nordic Lightning w picks, L 80 cm, left
500173	Hockey Stick Carbon Nordic Lightning w picks, L 82 cm, right
500174	Hockey Stick Carbon Nordic Lightning w picks, L 82 cm, left
500175	Hockey Stick Carbon Nordic Lightning w picks, L 84 cm, right
500176	Hockey Stick Carbon Nordic Lightning w picks, L 84 cm, left
500177	Hockey Stick Carbon Nordic Lightning w picks, L 86 cm, right
500178	Hockey Stick Carbon Nordic Lightning w picks, L 86 cm, left
500179	Hockey Stick composite Warrior w picks, L 70 cm, right
500180	Hockey Stick composite Warrior w picks, L 70 cm, left
500181	Hockey Stick composite Warrior w picks, L 72 cm, right
500182	Hockey Stick composite Warrior w picks, L 72 cm, left
500183	Hockey Stick composite Warrior w picks, L 74 cm, right
500184	Hockey Stick composite Warrior w picks, L 74 cm, left
500185	Hockey Stick composite Warrior w picks, L 76 cm, right
500186	Hockey Stick composite Warrior w picks, L 76 cm, left
500187	Hockey Stick composite Warrior w picks, L 78 cm, right
500188	Hockey Stick composite Warrior w picks, L 78 cm, left
500189	Hockey Stick composite Warrior w picks, L 80 cm, right
500190	Hockey Stick composite Warrior w picks, L 80 cm, left
500191	Hockey Stick composite Warrior w picks, L 82 cm, right
500192	Hockey Stick composite Warrior w picks, L 82 cm, left
500193	Hockey Stick composite Warrior w picks, L 84 cm, right
500194	Hockey Stick composite Warrior w picks, L 84 cm, left

See next page for more info about the accessories.

*Pictures shown may vary from actual products.*



Seat belt with jacking mechanism.

## Hockey Sticks

There are two different types of hockey sticks to Hockey Proff.

### Hockey Stick, Nordic Lightning Carbon, with picks:



These are handmade in Norway. Hockey sticks Nordic Lightning Carbon will provide a very efficient power transmission, and weighs about 40 grams less than a hockey stick in composite.

### Hockey Stick Warrior in composite, with picks:



Carbon provides Warrior hockey stick the rigidity that gives the best transmission. This is combined with glass fiber to increase strength and durability of the hockey sticks.



## Pushbar

If you need help getting around on the ice you can mount a pushbar on the hockey sledge. It is angle adjustable and easy to mount/dismount.



## Support Skates

If You've never played sledge hockey before, You can use support skates that are adjustable in width. These will make it easier for beginners and severely disabled to keep their balance without falling over. These can easily be removed when they are no longer needed.

*Pictures shown may vary from actual products.*

Hockey Proff is CE-marked and TÜV-tested

