

HandiSnow-5

2017



HandiSnow-5 is CE-marked and TÜV-tested

HandiSnow-5 is a Norwegian produced cross country sledge of high quality. With a frame made of aluminum, the sledge weighs less than 7 kg.



The sledge is suitable for both one-legged and double amputees and persons with reduced function in their legs.

Optimal spike position

The seat plate and knee support is adjustable, so that they can be angled in both front and back to obtain an ideal seat angle. The adjusted seat angle and seat height optimizes a good spike position, while the kneeling position ensures a very good weight distribution between the seat and knees. HandiSnow-5 has a backrest and 2 seat belts with jacking mechanism that keeps user firmly on the sleigh.

HandiSnow-5 has a functional handbrake. It is equipped with a brake strap between the brake lever and the seat so that it is easy to reach. In this way, you can safely brake without getting out of balance.

By mounting a higher back and wide waist belt, HandiSnow-5 can be used by people with bigger disabilities.

The skis are easy to take on and off, making the sledge easy to transport. In addition to ski bindings, the skis are equipped with solid fastening hooks to ensure good attachment of the sledge.

Standard equipment:

- Backrest H-15 cm
- Seat cushion
- Knee support, D-18 cm
- Handbrake
- Seat belts with jacking mechanism for hips and thighs
- Belt for kneebrace with fixlock buckle
- Belt for legs
- Ski with bindings and fastening hooks
- Ski poles with Pro Fit loops for better grip

Technical data:

Model	Article no.	Seat width	Seat depth	Seat height
HandiSnow-5 SW30	503217-30	30 cm	25 cm	47.5 to 50 cm
HandiSnow-5 SW35	503217-35	35 cm	25 cm	47.5 to 50 cm

Model	Back height	Ski length	Total weight	User weight	Seat Height front
HandiSnow-5 SW30	15 cm	182 cm	6,5 kg	110 kg	43 – 46 cm
HandiSnow-5 SW35	15 cm	182 cm	7 kg	110 kg	43 – 46 cm

Accessories:

Art.no.	Description
503218	Backrest, complete, SW35, BH20 cm
503218-20	Backrest, complete, SW30, BH20 cm
503219	Waist belt, wide
503292	Knee plate aluminum, W28, D22 cm
503344	Kneepads, W28, D22, H4 cm
503345	Kneepads, W28, D22, H2 cm
503347	Kneepads, W28, D18, H2 cm
500147	Poles with asphalt/ice spikes, L 120 cm, Pro Fit loops, size M, pair
500148	Poles with asphalt/ice spikes, L 125 cm, Pro Fit loops, size M, pair
500149	Poles with asphalt/ice spikes, L 130 cm, Pro Fit loops, size L, pair
500150	Poles with asphalt/ice spikes, L 135 cm, Pro Fit loops, size L, pair
500151	Poles with asphalt/ice spikes, L 140 cm, Pro Fit loops, size XL, pair
500152	Poles with asphalt/ice spikes, L 145 cm, Pro Fit loops, size XL, pair
503269	Seat belt M2 with jacking mechanism, Small
503270	Seat belt M2 with jacking mechanism, Medium
503271	Seat belt M2 with jacking mechanism, Large
503277	RollerSafe rollerskis with hydraulic brakes (poles not included)
503358	Cover for knees/thighs from the cold/rain
503310	Protective bag with straps



Seat belt with jacking mechanism.

Unique accessories

HandiSnow-5 can be used with RollerSafe rollerskis on permanent and regular basis during the summer. RollerSafe is designed and produced in Norway. These are the first rollerskis on the market with mounted hydraulic brakes operated via Bluetooth signals from brake levers on the poles. This gives the user full control and extra safety when braking and going downhill.



HandiSnow-5 with RollerSafe rollerskis.



RollerSafe rollerskis with hydraulic brakes (poles are not included). Delivered in a practical suitcase.

Pictures shown may vary from actual products

HandiSnow-5 is CE-marked and TÜV-tested

